

Safety Precautions

Congratulations on purchasing the Lux Flow Ionic Foot Bath! To ensure safe and enjoyable use, please follow these guidelines:

Before Use

1. Read this manual carefully and understand the foot bath's operation and safety features.
2. Inspect the foot bath for any damage or wear and tear before each use.
3. Ensure proper cleaning and maintenance according to the manual's instructions.

During Use

1. Remove any metal jewelry, especially from feet and ankles, to avoid electrical conductivity.
2. Avoid using the foot bath with open wounds, sores, or active infections.
3. Keep the foot bath away from children and pets.
4. Follow recommended usage guidelines and duration.
5. Handle the power cord and electrical components with care: avoid pulling, yanking, or stretching the wires, as this may cause damage or electrical shock.

Potential Hazards

1. Electrical shock: Avoid exposure to water or moisture beyond the designated water line.
2. Skin irritation: Discontinue use if you experience any discomfort, redness, or irritation.

Special Considerations

1. Pregnancy or breastfeeding: Consult your healthcare provider before use.
2. Pacemakers or implants: Consult your healthcare provider before use.
3. Sensitive skin: Start with shorter sessions and monitor your skin's response.

After Use

1. Clean and dry the foot bath according to the manual's instructions.
2. Store the foot bath in a safe and dry place, with the power cord properly wrapped and secured.